



vNAT[®] Buffer

- ✓ Inactivates viral pathogens in **5 minutes**
- ✓ Extracts and preserves viral nucleic acids
- ✓ Reduced assay time
- ✓ Reduced instrumentation cost (Cost Effective)



PRODUCT
vNAT[®] Buffer

PRODUCT SIZE
500 rxns

CATALOG NO
T00008

BioeXsen
Health Technology

BIOEXSEN GmbH
Brandshofer Deich 10, 20539, Hamburg, GERMANY
Phone: +49 40 228 58 56 0, Fax: +49 40 228 58 56 99
Web: www.bioexsen.com, e-mail: info@bioexsen.com